



Method for Coupling IgG to Microsphere Technology FERROSPHERE-N

Activation:

- 1) Weigh out 5mg of dry amine-functionalised FERROSPHERES-N (MTL/FER) into a 1.5mL microcentrifuge tube.
- 2) Wash 6 times (1mL/wash) with 10mM phosphate buffered saline, pH 7.2, 150mM sodium chloride (PBS), using a Dynal MPC™-S magnetic separator, or similar. Remove and replace the liquid carefully each time using a 1mL pipette.
- 3) Add 1mL of 10% glutaraldehyde in PBS (we recommend the use of Polysciences, Inc. 50% glutaraldehyde (Cat #: 18428), because of its high purity/low polymer content).
- 4) Incubate at room temperature for 2 hours with constant agitation (we recommend using a rotating blood tube wheel (Stuart rotator SB3, or similar) set to 30 rpm to ensure good, end over end, mixing).
- 5) Wash the FERROSPHERES-N 6 times with PBS as above.

Antibody Coupling:

- 6) FERROSPHERES-N have an IgG coupling capacity of 2-4mg/g. We recommend that IgG be added to the FERROSPHERES-N at a 5 to 10 fold excess of this coupling. Add the IgG to the FERROSPHERES-N in a total volume of 1mL, the extra volume being made up with PBS.
- 7) Incubate at room temperature for 18 hours with constant agitation on the tube rotator (for this extended incubation period, we recommend a setting of 15rpm).
- 8) Use the magnetic separator to pull the FERROSPHERES-N to the side of the tube and add 25 μ L of 1M sodium cyanoborohydride (dissolved in deionised water) to the FERROSPHERES-N/PBS mixture. Cyanoborohydride is a reducing agent and as such may impact on the activity of some antibodies. It is advisable to assess this beforehand.
- 9) Remove the tube from the magnetic separator and incubate at room temperature with constant agitation on the tube rotator (15rpm) for a further 30min.
- 10) Use the magnetic separator to pull the FERROSPHERES-N to the side of the tube and add 45 μ L of 100mg/mL adipic dihydrazide (adipic acid hydrazide), dissolved in deionised water, to the FERROSPHERES-N/PBS/cyanoborohydride mixture. The dihydrazide may require a little warming to dissolve in water at this concentration.
- 11) Remove the tube from the magnetic separator and incubate at room temperature with constant agitation on the tube rotator (15rpm) for a further 30min.
- 12) Wash the FERROSPHERES-N 6 times with PBS as above.

Removal of non-specifically bound antibody:

- 13) Wash the FERROSPHERES-N 3 times with deionised water.
- 14) Wash the FERROSPHERES-N 1 time with 50mM sodium acetate, pH 4.0, containing 300mM sodium chloride.
- 15) Wash the FERROSPHERES-N 1 time with deionised water.
- 16) Wash the FERROSPHERES-N 1 time with 50mM sodium carbonate, pH 9.0, containing 300mM sodium chloride.
- 17) Wash the FERROSPHERES-N 1 time with deionised water.
- 18) Repeat steps 14-17 a further 2 times.
- 19) Wash the FERROSPHERES-N with the storage buffer of choice and store in this buffer at 4 degrees Celsius.

Note: Capacity measurements were carried out by assaying the IgG coupled FERROSPHERES-N after the final water wash, using the Pierce BCA assay

The procedure described can be scaled up as required, but would require the use of a Dynal MPC™-1 magnetic separator (or similar).